	NON STERILE MEDICAL DEVICES Abdoma Abdominal gauze swab Absorba Packing Matocomp Gauze swabs Matovlies Nonwoven swabs Seton Gauze roll Tupfer Gauze tampon
DESCRIPTION Surgically invasive dressing materials made of bleached cotton gauze, threads and cotton loops, polyester threads, X-ray detectable thread and/or viscose-polyester nonwoven.	
INTENDED PURPOSE Devices: Absorba Packing, Matocomp Gauze swabs, Matovlies Nonwoven swabs, Seton Gauze roll, Tupfer Gauze tampon after sterilisation are designed to absorb liquids during surgeries and ambulatory procedures. Abdoma Abdominal gauze swabs after sterilisation are designed to absorb liquids and protect organs during surgical procedures.	
STERILITY The device is delivered non sterile. Recommended to be sterilised by hot steam in accordance with EN ISO 17665-2 or ethylene oxide in accordance with EN ISO 11135, after repacking under controlled conditions into sterilisation packaging, in accordance with EN ISO 11607-1 and EN ISO 11607-2.	
PREPARATION FOR USE 1. Before use, check that the package is not damaged. 2. After opening the package, check if the device is not dirty or damaged. 3. Subject the device to the sterilisation process as indicated. 4. Use in accordance with the intended use.	
CONTRAINDICATIONS Do not use the device for indications other than those specified in this instruction.	
PRECAUTIONS 1. Do not use directly on a wound without using a validated method of sterilisation. 2. Use the devices only in the form and size provided. 3. The device cannot be used after the expiry date. 4. Do not use it if the package is opened or damaged. 5. The device is intended for single use only. A reuse of the device is associated with a very high risk of patient infection and is unacceptable.	
STORAGE Storage and transport conditions are in accordance with the information in document "Summary of storage and transport conditions for TZMO SA finished products".	
UTILISATION Dispose of as medical waste in accordance with local regulations.	
REPORTING INCIDENTS Report medical incidents: qualitymanagement@tzmo-global.com	
DATE/LEAFLET VERSION Date of issue: 08_2024	

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