CZ

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ŽEBRA

EN

INSTRUKCIJS FOR USE
CODUBIX® S RIBS
Sterile prosthesis for filling of extensive defects in the chest wall

PRODUCT CHARACTERISTICS

- Codubix® S Ribs prostheses for filling of extensive defects in the chest wall are knitted of:
- polyester yarn – supplying sufficient strength against kinking and compression, as well as surface spread,
 - polypropylene yarn – featuring low specific weight and low melting point, which makes it possible to give the prosthesis adequate stiffness and hardness.
- High strength, low specific weight, hydrophobicity, contribute to the lack of absorption of liquids and physiological fluids; non-toxicity, chemical inactivity, high degree of healing in allow to identify the Codubix® S Ribs prostheses of extensive defects in the chest wall with the natural bone due to their physical properties. The Codubix® S Ribs prostheses can be easily modelled during the implantation and enable diagnostics without giving artefacts.

INDICATIONS

Codubix® S Ribs knitted polypropylene-polyester prostheses are designed to fill chest wall fragments after the treatment of primary and metastatic tumours affecting this structure and infiltrating non-small cell lung cancer.

CONTRAINDICATIONS

Do not use in people allergic to the ingredients of the product: polypropylene and polyester. It is not recommended to implant the Codubix® S Ribs in children and where the physiological growth process limits their use.

APPLICATION *)

Prior to Codubix® S Ribs prosthesis implantation the following must be specified:

- size of the defect,
- shape of the defect,
- curvature relevant to the patient's body.

Final match and adjustment of the prosthesis to the defect occurs intraoperatively (cutting of the prosthesis is possible). After it has been modelled the prosthesis needs to be carefully cleaned of the dust accumulated during the adjustment by rinsing in normal saline solution.

*) Developed by: Prof. Dr. med. M. Zawirski, Dr. med. P. Kasprzak, Department of Neurosurgery and Oncology of the Nervous System of the Medical University of Lodz. The cleaning procedure (time, quantity of saline solution) can vary depending on the shape and size of the cut prosthesis. Ensure thorough cleaning of all visible contaminants.

NOTE

- Codubix® S Ribs prostheses for filling of extensive defects in the chest wall are offered in a sterile version.
- The product should not be re-sterilised.
- Codubix® S Ribs must not be implanted after the expiry date.
- The product is intended for single use only. Due to the implantation method there is no possibility to reuse the product.
- Adjusting and cleaning of the prosthesis should be performed under aseptic conditions.
- Apogynic product.

POSSIBLE COMPLICATIONS

The implantation of bone implants can cause numerous types of complications i.e.:

- in the early post-operative period: infection, hematoma, swelling of the skin flap, development of blood and serous fluid exudates, increase in body temperature, ischaemia of the skin flap,
- in the long-term postoperative period: pain in the scar, sensation disorder, prosthetic collapse.

The aforementioned complications may contribute to the prolongation of treatment time, re-intervention of the surgeon, and in extreme cases, re-implantation of the prosthesis.

STORAGE CONDITIONS

Codubix® S Ribs prosthesis for filling of extensive defects in the chest wall should be stored:

- in dark rooms,
- at temperature range 15 + 35°C,
- at humidity 25 – 65%.

in the original unit packaging (box and direct packaging). Any mechanical damage to the packaging may expose the product to the loss of sterility.

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