
























SYMBOLS USED ON MEDICAL DEVICE LABELS

SYMBOL	MEANING
	CE with number – confirmation of meeting EU requirements
	manufacturer
	date of manufacture
	sterile product; sterilised using steam
	sterile product; sterilised using ethylene oxide
	sterile product; sterilised using irradiation
	sterile product; no information about the type of sterilisation when more than one method can be used
	non sterile product
	batch code – the manufacturer's code
	use-by-date – indicates the date after which the medical device is not to be used
	do not re-use – product for single use only
	do not re-sterilise
	do not use if package is damaged
	not allowed for food contact
	allowed for food contact
	keep away from sunlight
	keep dry
	temperature limits to which the medical device can be safely exposed
	product contains natural rubber latex
	product is latex-free
	product gives protection against microorganisms
	product gives limited protection against chemicals
	instructions for use attached to the product

SYMBOLS USED ON MEDICAL DEVICE LABELS



range of humidity to which the medical device can be safely exposed



„Green dot” means that the manufacturer of the product contributes to the cost of recovery and recycling



information about the raw material the given packaging is made of



throw away your trash



non-pyrogenic product



hand wash



wash at the indicated temperature



drip dry



do not iron



cannot be bleached with chlorine



do not dry clean



wash at normal temperature in the tumble drier



must not be bleached



not for contact with food



„hands” show from which side should the packaging be opened




this side up



handle with care

KJ

quality control approval



registration number for Ukraine - informs since when the product is registered



registration number for Russia - informs since when the product is registered

Per. № ФСЗ 2012/12667
от 22.08.2012 г.