

Made of transparent polyurethane film secured with a silicone paper liner.

The film edges have colourful “fingerlifts” enabling the application of the product to the designated operational area. The film has an additional silicone paper wrapping.

Characteristics:


- is flexible and extensible
- does not irritate the skin
- is latex-free
- does not cause allergic reactions
- is not cytotoxic
- is permeable to water vapour

Application:

- to protect the surgical field against microbial penetration

Size and packing:

- sterile (sterilised with ethylene oxide): packed in a box by 10 pcs
- wide range of sizes: from 20 cm x 10 cm to 90 cm x 60 cm



**MD**

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>