

- made of chlorine-free bleached cotton gauze
- edges are secured to prevent the release of threads
- have X-ray detectable thread, which is permanently placed inside the swab
- may have an additional element - cotton loop, which allows fixing the swab outside the OP-field
- can be sterilised

Application:

- for draining large amounts of body fluids during surgical procedures

Types:

- sterile and non sterile
- with X-ray detectable thread and additionally with or without a cotton loop
- number of layers: 2-, 4-, 6-, and 8-ply

Size and packing:

- sterile: packed by 1 pc, 2 pcs, 3 pcs, 5 pcs or 10 pcs
- non sterile: packed by 10 pcs, 20 pcs, 30 pcs, 45 pcs or 60 pcs
- available sizes: 30 cm x 30 cm, 30 cm x 45 cm, 45 cm x 45 cm, 45 cm x 70 cm, 50 cm x 50 cm, 75 cm x 90 cm



Devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>