

## Managing the infection

- vapour-permeable polyurethane film constituting the outer layer of the dressing allows for gaseous exchange and protects against external factors
- polyurethane foam layer absorbs excess exudate and protects the wound against mechanical damage
- active layer consisting of chitosan and silver provides an aseptic environment in the wound and creates favourable conditions for wound healing

## Application:

- moderate to high exuding wounds
- infected wounds or at risk of secondary infection
- shin ulcers
- pressure ulcer wounds
- postoperative wounds
- traumatic extensive injuries
- burn wounds
- wounds accompanying the so-called diabetic foot

## Types:

- sterile

## Sizes and packaging:

- 10 cm x 10 cm, packed by 5 pcs
- 10 cm x 20 cm, packed by 5 pcs



Sterile devices are certified by the Notified Body PCBC SA and bear the CE mark with number 1434

>>