

## Managing the infection

- the absorbent pad contains activated metallic silver which has antibacterial properties
- the dressing absorbs excess secretion and maintains a moist wound environment
- the elastic nonwoven fabric allows the product to adhere to the contours of the body shape
- steam permeable, breathable - allows for gas exchange

## Application:


- hard-to-heal wounds
- infected wounds and clean wounds as a dressing to prevent infections
- graft donor sites

## Types:

- sterile

## Sizes and packaging:

- 10 cm x 10 cm, packed by 5 pcs
- 10 cm x 20 cm, packed by 5 pcs



Sterile devices are certified by the Notified Body  
Superiore di Sanità and bear the CE mark with number  
0373

>>