

- made of cotton gauze and medical nonwoven, being the absorbent part
- the combination of gauze and nonwoven increases the absorption
- can be sterilised

Application:

- to absorb wound fluids during operations and ambulatory procedures


Types:

- sterile and non sterile

Size and packing:

- available sizes: 10 cm x 20 cm (folded), 9 cm x 9 cm (sewn)
- 10 cm x 20 cm sterile: packed by 1 pc, 2 pcs, 4 pcs
- 9 cm x 9 cm sterile: packed by 1 pc
- 10 cm x 20 cm non sterile: packed by 50 pcs

NOTE: For invasive procedures use the product only in an unwound state.



Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

