

Applicators are non sterile reusable surgical instruments used in the surgical treatment of urinary incontinence in women. Made of stainless martensitic steel meet the requirements of standards concerning surgical tools.

Application:

- as accessories used in surgical treatment of stress incontinence in women taking into account the division into methods applied – each tool is used to create a tunnel in patient’s body in order to properly implant the urological tape

Types:

- TVT applicator is an accessory used in TVT method (Tension-free Vaginal Tape)
- TOT applicator is an accessory in TOT method (Trans Obturator Tape)
- a set of applicators for TOT consists of two applicators which are mirror images. It is used in TOT method (Trans Obturator Tape)
- a set of applicators for TOT consists of two applicators which are mirror images. It is used in “outside-in” TOT method (Trans Obturator Tape)
- non sterile (steam sterilisation is recommended)

Packing:

- TOT and TVT applicators packed by 1 pc
- set of TOT applicators packed by 2 pcs

More information on: www.dallop-nm.com



MD

>>