

- components of plastic (e.g. [umbilical clamps](#), forceps etc.) or metal in sterile packaging
- single use symbol engraved on all metal instruments
- sharp instruments are secured with a special cap, there is no risk of puncturing the packaging
- sterilised in validated process (in accordance with EN ISO 17665-1 and EN ISO 11135)
- instruments cannot be re-sterilised

Application:

- for use during surgical and outpatient procedures

Types:

- sterile
- forceps, scissors, tweezers, scalpel holders, etc.

Size and packing:

- sterilisation packaging by 1 pc; in addition, individual sterile instruments can be packed in a cardboard dispenser enabling the organisation and storage of products in a hospital



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>