

- 100% cotton
- chlorine-free peroxide bleached
- no loose fibres on the gauze surface
- high capacity to absorb moisture and liquids
- can be sterilised

Application:

- for general wound care (sterile)
- for hygienic and household purposes (non sterile)

Types:

- sterile (steam sterilised) and non sterile
- 13-, 17-thread

Size and packing:

- sterile in sachets: 1 m²; 1/2 m²; 1/4 m²
- non sterile in sachets: 1 m²; 1/2 m²
- folded: width from 0.90 m to 1.20 m, lengths: 50 m, 100 m, 200 m

NOTE: For invasive procedures use the product only in an unwound state.



Sterile devices are certified by the Notified Body TÜV Nord Polska sp. z o.o. and bear the CE mark with number 2274

