

Cohesive fixing bandage with min. 85% elasticity. The bandage is self-adhesive so there is no need to use any clips for fixing. It conforms perfectly to the body shape, and gives a reliable and non-slip fit. The bandage adheres only to itself, but not to skin, hair or clothing. It can be torn into pieces without the need to use the entire bandage at once.

Disposable.

Contains synthetic latex.

Application:

- for fixing dressings, especially the rounded body parts

Types:

- non sterile

Size and packing:

- non sterile: 8 cm, 10 cm and a length of 4.5 m – length measured in a stretched condition



MD

>>