

Cohesive fixing bandage with 90% elasticity.

Latex-free!

The bandage is self-adhesive so there is no need to use any clips for fixing.

Conforms perfectly to the body shape, and gives a reliable and non-slip fit. The bandage adheres only to itself, but not to skin, hair or clothing.

Disposable.

Application:

- for fixing dressings, especially the rounded body parts

Types:

- non sterile
- sterile

Size and packing:

- non sterile: 6 cm, 8 cm, 10 cm, 12 cm and a length of 4 m – length measured in a stretched condition



**MD**

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>