

Used to fix dressings on the body in lieu of a traditional bandage, especially on body parts difficult to dress, such as the pelvis or groin.

They are in form of a sleeve, making it easy to apply - there is no need to know the bandaging techniques to put on the mesh.

8 different sizes - every size of the circuit is adjusted to a specific body parts.

If chosen appropriately the elastic structure will make the mesh fit perfectly.

Quick to apply, which is especially important when changing dressings in children.

Ideal for home use - keep in mind that inappropriately applied bandages can disturb circulation - a single-layer construction and 8 sizes to choose from minimise that risk.

Application:

- to fix dressings on difficult-to-dress areas, such as: pelvis, perineum
- for dressing surgical wounds of the abdomen, chest and head
- for supporting dressings on burn wounds
- in paediatric surgery

Types:

- non sterile

Size and packing:

- 1 finger
- 2 hand, finger
- 3 hand, foot
- 4 lower leg, knee, arm, foot, elbow
- 6 head, arm, lower leg, knee
- 8 thigh, head, hips
- 10 hips, abdomen
- 14 chest, abdomen, hips
- available lengths: 1 m and 10 m (respectively: 2.5 m, 25 m when stretched)
- packed by 1 pc



MD

