

- made of porous transparent air and water vapour permeable film
- allows observation of the fluid level in fixed probes, cannulas and catheters
- very gentle to the skin; can be removed painlessly without leaving any residues
- does not absorb moisture
- serrated edges on both sides for easier tearing off the desired length in the hand
- acrylic adhesive does not cause skin irritation

## Application:

- for supporting dressings of the face and other sensitive areas
- for fixing tubes, catheters and other medical devices

## Types:

- non sterile

## Size and packing:

- available lengths: 5 m and 9.14 m and a width of 1.25 cm, 2.5 cm and 5 cm
- non sterile: packaging depending on the size



**MD**

>>