

Made of transparent polyethylene film.


The elasticity and softness of the material allows for easy and safe manipulation during surgery. It is equipped with an adhesive tape and a cardboard box with cut-outs, so that the cover can be easily fixed on the device.

Application:

- to cover medical devices, universal for all cameras and endoscopes, fibre optic and cables

Size and packing:

- sizes: 13 cm x 235 cm, 13 cm x 244 cm, 15 cm x 250 cm
- sterile (sterilised with radiation or ethylene oxide): packed in a paper-film by 1 pc



**MD**

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>