

Made of polyethylene or polypropylene film.

The covers (except for the lamp handle covers) have an attached elastic band or adhesive tape facilitating application.

Application:

- to drape medical devices

Types:

- for equipment
- for C-arm
- for microscope
- for lamp handle

Size and packing:

- wide range of sizes and shapes
- sterile (sterilised with ethylene oxide): packed in a paper-film pouch by 1 pc



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>