

Special packaging for radiation sterilisation.

Application:

- as packaging for sterilisation by radiation (R) and ethylene oxide (EO) with an additional indicator for the sterilisation method

The product meets the requirements of EN 868-6, EN 868-5, EN ISO 11607-1 and EN ISO 11607-2.

Types:

- pouches
- reels

Size and packing:

- pouches - range of sizes: width from 40 mm to 600 mm; length from 90 mm to 1200 mm; packed by 100 pcs
- reels - range of sizes: width from 50 mm to 650 mm; length of 60 m and 200 m; packed by 1 pc



MD

>>