

- the outer layer is made of chlorine-free bleached cotton gauze
- the core is made of absorbent cotton wool
- the oval form reduces its size and facilitates the use in body cavities
- the attached cotton thread makes it easy to remove the tampon from the body
- can be sterilised

Application:

- for gynaecological procedures

Types:

- sterile

Size and packing:

- available sizes: 30 cm x 30 cm
- sterile: packed by 1 pc

NOTE: For invasive procedures use the product only in an unwound state.



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>