

- made of 100% cotton gauze, chlorine-free peroxide bleached
- no loose fibres on the surface
- edges secured to prevent release of any fibres
- high capacity to absorb moisture and liquids
- can be sterilised

Application:

- to absorb fluids during outpatient procedures and operations (sterile)
- for the hygienic and cosmetic purposes (non sterile)

Types:

- sterile (steam sterilised) and non sterile
- 13-, 17-thread
- with or without X-ray detectable thread
- tied in piles or loosely arranged
- 8-, 12-, 16-ply

Size and packing:

- non sterile: 5 cm x 5 cm, 7.5 cm x 7.5 cm, 10 cm x 10 cm, 10 cm x 20 cm packed in a paper bag by 100 pcs
- sterile: 5 cm x 5 cm, 7.5 cm x 7.5 cm, 10 cm x 10 cm, 10 cm x 20 cm packed in blister packaging by 2 pcs, 3 pcs, 5 pcs

NOTE: For invasive procedures use the product only in an unwound state.



MD

Devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>