

- the outer layer is made of medical nonwoven
- the inner layer is made of cellulose pulp
- the use of cellulose pulp increases the absorption property, making it the most absorption swab
- can be sterilised

## Application:

- for dressing wounds and heavily bleeding injuries

## Types:

- sterile and non sterile

## Size and packing:

- available sizes: 10 cm x 15 cm, 10 cm x 20 cm
- sterile: packed by 1 pc
- non sterile: packed by 20 pcs

NOTE: For invasive procedures use the product only in an unwound state.



Sterile Devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

