

Table cover with 2 divided areas meets the requirements of standard EN 13795.

The absorbent layer in the middle part of the cover collects liquids, and the film layer stops them and prevents them from soaking through onto the table.

It is sold separately, as a set element or as an integral part of OP-sets (set wrapping).

Application:

- to cover an instrument table, on which drapes or instruments are placed

Size and packing:

- size: 150 cm x 90 cm
- sterile (sterilised with ethylene oxide): packed in a paper-film pouch by 1 pc; in addition, individual sterile covers can be packed in a cardboard dispenser enabling the organisation and storage of products in a hospital



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>