

Combination of comfort and safety.

Latex examination gloves are characterised by high elasticity and mechanical strength. Provide an effective protective barrier, and a beaded cuff enables efficient and secure donning.

Ambulex – examination and disposable protective, powdered

- powdered with natural corn starch, which insulates the hand from latex layer
- textured fingertips provides a better grip for efficient instrument handling
- ambidextrous – fit both hands

Ambulex gloves classified as personal protective equipment (PPE) of category III, comply with EN ISO 374. Before using the gloves refer to the [add. information for users](#). For detailed information on chemical resistance, please contact your medical sales representative.

Ambulex P – examination and disposable protective gloves, powder-free – especially recommended for people with sensitive skin

- made from latex, powder-free
- polymer coating protects enables effortless damp and dry hand donning
- rough surface provides a secure grip for efficient instrument handling in dry and wet conditions
- ambidextrous – fit both hands

Ambulex P gloves classified as personal protective equipment (PPE) of category III, comply with EN ISO 374. Before using the gloves refer to the [add. information for users](#). For detailed information on chemical resistance, please contact your medical sales representative.

[Ambulex P – Declaration of Conformity](#)

Application:

- widespread use in healthcare
- reduce the risk of infection during the study and treatment, and contact with contaminated medical materials

Types:

- non sterile
- powdered and powder-free

Size and packing:

- XS, S, M, L, XL in a dispenser box by 100 pcs
- powdered Ambulex gloves size M packed in a box by 10 pcs dedicated for household use



MD

>>