

Practical gynecological workshops 29th March 2023

Lightweight urological tapes are the latest generation of sling implants used in the treatment of female stress urinary incontinence. They are knitted of transparent and blue monofilament polypropylene yarn. Blue orienting lines along the tape facilitate the implantation process and increase the visibility of the implant in the operating field.

On both endings the tape is equipped with blue loops which allow safe passing through anatomical structures of the patient. The handles are also made of monofilament polypropylene yarn.

Due to the ultralight surface mass, atraumatic edges, excellent strength parameters, high elasticity and very good biocompatibility, the product can be safely and effectively used in the surgical treatment of urinary incontinence.

Application:

- urological tapes are recommended for use during the surgical treatment of Stress Urinary Incontinence (SUI). Surgical treatment of urinary incontinence involves suspending the urethra with a tension-free urological tape.
- the device is intended for implantation by the following methods: Tension Free Vaginal Tape (TVT), Trans Obturator Tape (TOT) or Tension Free Vaginal Tape-Obturator (TVT-O)



Types:

- sterile (sterilised with ethylene oxide)

Size and packing:

- length of 30 cm and 45 cm
- packed in double paper-film packaging with gusset and then in a cardboard box by 1 pc

More information on: www.dallop-nm.com



MD

Sterile devices are certified by the Notified Body PCBC SA and bear the CE mark with number 1434

>>