

The nonwoven layer on the top allows to absorb fluids from instruments used during an operation, and the film part provides barrier to liquid penetration.

The telescopic folding facilitates the application on the table and protects against accidental contamination while the sandblasted film prevents the cover from slipping.

Complies with EN 13795.

Application:

- to cover Mayo tables

Types:

- standard
- with thicker film

Size and packing:

- size: 145 cm x 80 cm
- sterile (sterilised with ethylene oxide): packed in a paper-film pouch by 1 pc; in addition, individual sterile covers can be packed in a cardboard dispenser enabling the organisation and storage of products in a hospital



Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

