

The elasticity of the raw material, makes it easy to adjust to any shape, and the applied “fingerlifts” allow for quick application.

The Velcro tape is composed of two parts: hooks and loops. Both parts have a glue layer and a protective release paper with “fingerlifts”. One part of the tape is to be attached to a given place on a drape, and the second part can be removed or attached repeatedly again as often as needed.

Application:


- for fixing drapes, tubes, cables

Types:

- with or without Velcro

Size and packing:

- sterile (sterilised with ethylene oxide): packed in a paper-film pouch by 1 or 2 pcs
- sizes: 50 cm x 9 cm, 50 cm x 10 cm, 90 cm x 5 cm, 75 cm x 10 cm (without Velcro) and 25 cm x 2 cm (with Velcro)



**MD**

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>