

Medical surgical caps are made of thin polypropylene nonwoven. They can have ties for fixing or a delicate elastic band that holds the cap on the head. The offer also includes medical caps with a sweatband.

Application:

- in operating theaters and outpatient clinics to protect against transmission of infectious agents between the patient and medical staff during surgical procedures

Types:

- with elastic band, with ties, with sweatband
- sterile (steam sterilised) and non sterile
- colour: blue

Size and packing:

- sterile: packed in blister packaging by 1 pc
- non sterile: packed in a cardboard box by 100 pcs
- non sterile: packed in bag of 100 pieces (applies to a cap with ties and a sweatband)



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>