

Thinking of people who are allergic to latex.

Nitrile examination gloves are characterised by the absence of natural rubber latex and high tensile strength. Provide an effective protective barrier, and a beaded cuff enables efficient and secure donning.

Ambulex NITRYL - nitrile examination gloves, powder-free

- powder-free
- high tensile strength both at the beginning and end of the validity period
- textured fingertips provide a secure grip for efficient instrument handling in dry and wet conditions
- ambidextrous - fit both hands
- available colours: white, violet, black (especially for tattoo studios), elongated white (an elongated cuff allows for better protection against infections)

The powder-free nitrile gloves are classified as personal protective equipment of category III, comply with EN ISO 374. Before using the gloves refer to the [add. information for users](#). For detailed information on chemical resistance, please contact your medical sales representative.

Application:

- widespread use in healthcare
- reduce the risk of infection during examinations and procedures, and contact with contaminated medical materials

Types:

- non sterile
- powder-free

Size and packing:

- white, violet, elongated white: XS, S, M, L, XL
- black: S, M, L
- in a dispenser box by 100 pcs



MD

>>