

- made of white nonwoven fabric which easily adapts to the shape of the body
- porous structure is extremely permeable to both air and water vapour, thereby allowing the skin to breathe and facilitating evaporation of perspiration
- serrated edges on both sides for easier tearing off the desired length in the hand
- acrylic adhesive does not cause skin irritation, enhances permeability to air and is particularly well tolerated by the skin

Application:

- for securing a dressing on a sensitive part of the body
- for fixing probes, catheters and other medical devices in people with hypersensitive skin

Types:

- non sterile

Size and packing:

- available lengths: 5 m and 9.14 m and a width of 1.25 cm, 2.5 cm and 5 cm
- non sterile: packaging depending on the size



MD

>>