

- backing material made of white elastic nonwoven, thanks to which the dressing conforms well to the cannula and the patient's body shape
- porous nonwoven allows the skin to breathe
- special notch for the cannula for easy application
- micro-net prevents the dressing from sticking to the wound
- additional absorbent pad under the foot of the cannula
- applied adhesive does not cause skin irritation
- rounded edges protect the dressing against premature detachment from the patient's body
- wide protective paper overlap for easy dressing application even when wearing gloves

Application:

- for the sterile wound treatment of the puncture point and securing the indwelling cannula at the same time

Types:

- sterile

Size and packing:

- sterile: packed by 1 piece in a box by 50 pcs
- sizes: 5.8 cm x 8 cm, 7.2 cm x 5 cm



MD

Sterile devices are certified by the Notified Body
Superiore di Sanità and bear the CE mark with number
0373

>>