

FIXOPORE S nonwoven dressing with absorbend pad, FIXOPORE S nonwoven eye dressing with absorbend pad

- backing material made of white elastic nonwoven, thanks to which the dressing conforms well to the body shape and does not constrain moves
- porous nonwoven allows the skin to breathe
- absorbent core is located centrally and the surrounding adhesive holds the dressing firmly in place
- micro-net prevents the dressing from sticking to the wound
- applied adhesive does not cause skin irritation
- rounded edges protect the dressing against premature detachment from the patient's body
- wide protective paper overlap for easy dressing application even when wearing gloves
- X-ray permeable dressing

Application:

- the treatment of wounds and other surgical wounds, cuts and abrasions

Types:

- sterile
- shapes: rectangular, oval (after eye surgery)

Size and packing:

- sterile: packed by 25 pcs, 50 pcs or 100 pcs depending on the size
- rectangular available in sizes: 7.2 cm x 5 cm, 10 cm x 6 cm, 10 cm x 8 cm, 15 cm x 8 cm, 20 cm x 10 cm, 25 cm x 10 cm, 30 cm x 10 cm, 35 cm x 10 cm
- oval available in size: 6.5 cm x 9.5 cm



MD

Sterile devices are certified by the Notified Body
Superiore di Sanità and bear the CE mark with number
0373

>>