

- made of swab nonwoven fabric
- no loose fibres on the surface
- compact structure allows cutting without releasing any loose thread
- high capacity to absorb moisture and liquids
- very soft, perfectly fits to the body shape
- can be sterilised

Application:

- to absorb fluids during outpatient procedures and operations (sterile)
- for the hygienic and cosmetic purposes (non sterile)

Types:

- sterile (steam sterilised) and non sterile
- 30 g/m², 40 g/m²
- with or without X-ray detectable thread
- tied in piles or loosely arranged
- 4-, 6-, 8-ply

Size and packing:

- non sterile: 5 cm x 5 cm, 7.5 cm x 7.5 cm, 10 cm x 10 cm, 10 cm x 20 cm packed in a paper bag by 100 pcs
- sterile: 5 cm x 5 cm, 7.5 cm x 7.5 cm, 10 cm x 10 cm, 10 cm x 20 cm packed in blister packaging by 2 pcs, 3 pcs, 5 pcs

A special type of nonwoven swabs are tracheotomy swabs:

- made of swab nonwoven
- with O- or Y-shaped cut-out
- very soft, perfect for using with tubes, drains, catheters

Application:

- to secure tracheotomy tubes, catheters and other types of insertions

Size and packing:

- sterile: 5 cm x 5 cm, 7.5 cm x 7.5 cm, 10 cm x 10 cm packed in blister packaging by 2 pcs



Devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

