

Made of tear resistant SMS nonwoven in form of a tunnel. The covers are fixed on a device by means of ties or an elastic band.

Application:

- to protect hospital equipment and cables of medical devices against damage and contamination

Types:

- C-arm tunnel cover
- cable tunnel cover
- with ties or elastic band

Size and packing:

- sterile (steam sterilised): packed in a paper-film pouch by 1 pc
- non sterile: packed in a plastic bag by 70 or 100 pcs depending on the size
- sizes: 120 cm x 80 cm, 150 cm x 10 cm, 200 cm x 26 cm, 200 cm x 32 cm, 250 cm x 15 cm



**MD**

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>