

- the inside is made of swab nonwoven
- the swabs are secured with a knitted bandage, whose ends form ties
- the ties are tied at the head's back to fix the dressing under the nose

## Application:

- to absorb fluids after laryngological procedures

## Types:

- unsterile

## Size and packing:

- available sizes: 5 cm x 2.5 cm, 7.5 cm x 4 cm
- unsterile: packed by 50 pcs



**MD**

Sterile devices are certified by the Notified Body TÜV  
Nord Polska Sp. z o.o. and bear the CE mark with  
number 2274

>>