

Nursing surgical caps are made of thin polypropylene nonwoven. The fitted cut of the caps allows you to cover the head perfectly, and the flexible and delicate elastic band holds the cap on the head.

Application:

- in operating theaters and outpatient clinics to protect against transmission of infectious agents between the patient and medical staff during surgical procedures

Types:

- bouffant cap
- sterile (steam sterilised) and non sterile
- colour: blue

Size and packing:

- sterile: packed in blister packaging by 1 pc
- non sterile: packed in a cardboard box by 100 pcs



**MD**

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>