

Patient gowns are made of soft polypropylene nonwoven. They have short sleeves, and the neckline is finished with a pleat, which is also used to tie the product. The front ties are sewn into the gown to be tied around the waist. Our offer includes gowns for adults as well as gowns for children in small sizes.

Application:

- as patient's clothing in treatment rooms

Types:

- non sterile
- colour: navy blue

Size and packing:

- size: universal, 54/98, 104/128, 132/146
- non sterile: packed in a plastic bag by 10 pcs
- [dimensions of the final product](#)



**MD**

>>