

Made of nonwoven or laminate that meet the requirements of standard EN 13795. The covers are fixed on the limb by means of medical tapes.

Application:

- to cover legs during operations

Types:

- packed individually or are a part of OP-sets
- with or without added medical fixing tape
- raw material: SMS, TF, Blue Special and Blue Comfort

Size and packing:

- sterile (sterilised with ethylene oxide or steam): packed in a paper-film pouch by 1 pc
- sizes: from 60 cm x 35 cm to 130 cm x 140 cm



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>