

Control of highly exuding wounds.

- the absorbent properties of the dressing allow it to stay for a long time on wounds with a large amount of exudate and regulate the humidity by maintaining a moist environment that creates optimal conditions for the treatment process
- protection against mechanical injury - a layer of soft and flexible polyurethane foam minimises the effects of accidental secondary injury to a wound

Application on many types of exuding wounds, especially:

- burn wounds
- postoperative wounds
- chronic venous leg ulceration
- extensive traumatic injuries
- bedsores

Types:

- sterile

Size and packing:

- sterile: 10 cm x 10 cm, packed by 5 pcs



Sterile devices are certified by the Notified Body PCBC SA and bear the CE mark with number 1434

>>