

Accelerated healing and infection control

- resorbable dressings made of a modified polymer of natural origin – chitin, obtained from the shells of marine crustaceans
- thanks to the resorption capacity of the dressing, it does not need to be changed
- the breakdown products of the dressing in the wound together with proteins form a building material source for new cells
- regulates wound moisture and maintains a moist environment
- the silver content has an antibacterial effect (when using the Medisorb R Ag dressing)

Application:

It is intended for use in the case of chronic exudative wounds that are difficult to heal, infected and not infected, including:


- wounds due to venous insufficiency
- diabetic foot
- bedsores
- traumatic wounds
- surgical wounds
- burn wounds
- the POWDER version allows it to be used on deep wounds

Types:

- sterile
- Medisorb R (version MEMBRANE and POWDER)
- Medisorb R Ag

Sizes, versions and packing:

- Medisorb R: version MEMBRANE ø 9 cm; packed by 5 pcs, version POWDER: 1 g; packed by 5 pcs
- Medisorb R Ag: ø 9 cm; packed by 5 pcs



MD

>>