

A range of flat paper-film packaging for sterilisation.

High-quality materials used in the manufacture of the packaging as well as validated technological manufacturing process provide excellent “peel” effect – dust-free opening of the packages. This minimises the risk of contamination by cellulose particles, which can be emitted while opening paper packaging.

Due to the adhesive layer the pouch can be closed without the need to use a heat-sealer.

Chemical indicators of the sterilisation process are printed on the paper within the seal.

The product meets the requirements of EN 868-3, EN 868-5, EN ISO 11607-1, EN ISO 11607-2 and EN ISO 11140-1.

[Instruction for closing the pouches](#)

Application:

- as packaging for sterilisation by steam (STEAM) and ethylene oxide (EO)

Types:

- pouches

Size and packing:

- range of sizes: width from 40 mm do 60 mm; length from 90 mm to 1200 mm; packed by 100 pcs



**MD**

>>