

As part of disposable bedclothes, we offer:

- Vlieskomfort Set of nonwoven bedclothes
- Vlieskomfort Nonwoven pillow case
- Vlieskomfort Nonwoven duvet cover
- Vlieskomfort Nonwoven sheet and film coated nonwoven sheet

Disposable bedclothes are made of nonwoven polypropylene. Bedclothes are offered both as bedclothes sets, as well as individual components. Pillow cases and duvet covers have flaps that prevent the duvet, pillow or blanket from sliding out of the inside. The sheets may have elastic band sewn in for better protection of the mattress. The bedding is turned inside out to make it easy and quick to put the cover and pillow case on the duvet or pillow.

Application:

- as disposable bed linen in hospitals and surgeries

Types:

- non sterile and sterile
- colour: blue

Size and packing:

- [sizes and components](#)
- non sterile:
  - set of nonwoven bedclothes packed in plastic bag by 1 set
  - nonwoven pillow cases and duvet covers packed in plastic bag by 5 pcs
  - nonwoven sheets packed in plastic bag by 10 pcs
- sterile: packed by 1 set in paper-film pouch



Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>