

The shape and size of drape allow to drape the patient ensuring the safety for him/her and the medical staff during various procedures. The adjusted construction of the drapes makes it possible to quickly and effectively drape the patient and set the boundaries of the operating field. For even easier work the drapes are marked with pictograms indicating the unfolding directions.

The surgical drapes are made of the highest quality raw materials that meet the requirements of EN 13795, thus ensuring safety while using.

Some examples:


- angiography drape
- caesarean section drape
- hand/foot/limb drape
- gynaecological drape
- ophthalmology drape
- under buttocks drape for delivery
- urology drape

Types:

- drapes can be made of laminates (TF, Blue Basic, Blue Special, Blue Comfort) or nonwoven (SMS)
- to best adjust the drapes to a given operation, they can be equipped with additional elements, including adhesive tapes, openings, cut-outs, liquid collection bags, pouches, elastic elements, incise film, integrated legs, film panels, additional absorbent layers (so-called absorbent patches), cable organisers

Size and packing:

- wide range of sizes
- sterile, packed in sterilisation packaging by 1 pc; in addition, individual sterile drapes can be packed in a cardboard dispenser enabling the organisation and storage of products in a hospital



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>