

Medidress surgical sets are made of fabrics characterised by softness and delicacy. A special weave of yarn used in the fabric gives high comfort in using the clothing. The product may have a bacteriostatic finish that prevents the multiplication of bacteria, and the carbon thread inserted into the fabric prevents the accumulation of electric charges. The sets consist of a sweatshirt and trousers, or a sweatshirt and a skirt. The sweatshirt has a yoke on the front, short sleeves and three sewn pockets, which can be decorated with coloured vertical trimmings at the customer's request. On the sides of the sweatshirt there are slits that do not restrict the user's movements. Straight pants, finished with a belt with ties to adjust the waist. Straight skirt (envelope cut), with a belt and ties to tie at the waist.

Application:

- as clothing for medical staff in the operating unit, the product complies with EN 13795

Types:

- non sterile
- of polyester textile: blue, green, orange, grey
- of cotton-polyester textile: blue, green, red

Size and packing:

- size: XS, S, M, L, XL, XXL, 3XL; height: 158-164; 170-176; 182-188 - identify the size with the colour of the hanger in the top and of the straps in the trousers and the skirt
- packed by 1 set



MD

>>