

The Medidress Perfect surgical gowns are made of two types of fabrics: a light polyester fabric and a three-layer laminate in the critical zones of the product. The design of the gown allows it to be put on aseptically, and the long polyester cuffs fit tightly around the wrist. In addition, the cuff has a sewn-in loop to put on the finger, which prevents the sleeves from pulling up, which guarantees that the sleeve does not slip out from under the glove.

Application:

- as clothing for medical staff, used for long and highly invasive medical procedures, the product complies with EN 13795

Types:

- non sterile
- sterile - in the service offered by CitoNet companies
- colour: green, blue

Size and packing:

- S, M, L, XL, XXL
- packed by 1 pc



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274. The conformity assessment of products was limited to brand new products (it does not include the sterilization process after being placed on the market)

