

The gowns are made of polypropylene nonwoven and in the critical zones (front and forearms) they are reinforced with barrier laminate. Provide complete resistance to microbial and liquid penetration. Every gown is loose, wide, and enables to precisely cover the surgeon's trunk. The sleeves have a seamless dust-free polyester cuff which precisely embraces the wrist. Two sets of ties enable first "dirty" and then "clean" tying. Special cardboard with affixed external ties allows for aseptic application of the gown, and the book-folding enables quick and self-application by the surgeon.

Application:

- for long operations with a big amount of fluids

Types:

- sterile (sterilised with ethylene oxide) and non sterile
- sterile with nonwoven hand towels
- sterile with nonwoven hand towels and additional wrap in crepe paper
- colour: blue

Size and packing:

- sizes: M, L, XL, XXL
- sterile: packed in a paper-film pouch by 1 pc
- non sterile: packed in a plastic bag by 10 pcs



Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

