

- composition of sterile swabs of various sizes and number of threads
- whole package is adjusted to individual customer needs
- combination of swabs in a single package reduces the amount of waste
- one package eliminates the cost of the packaging of individual components
- one package eliminates the need to open multiple single in-between packaging, and thus does not increase the risk of infecting a sterile product

Application:

- to absorb fluids during surgical and outpatient procedures

Types:

- sterile
- gauze swabs: number of threads (13-, 17-thread), stratification (8-, 12-, 16-ply), with or without X-ray detectable thread
- nonwoven swabs: basis weight (30 g/m², 40 g/m²), stratification (4-, 6-, 8-ply), with or without X-ray detectable thread
- on request - tied up in piles for easy counting

Size and packing:

- combination of sizes: 5 cm x 5 cm, 7.5 cm x 7.5 cm, 10 cm x 10 cm, 10 cm x 20 cm contained in one package



MD

Sterile devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>