

- made of chlorine-free bleached cotton gauze with attached drawstrings
- the strings form the shape of the tampon and facilitate its removal from the OP-field
- can be sterilised

Application:

- for use during and outpatient procedures to control bleeding and absorb large amounts of fluids

Types:

- sterile
- with 1 or 2 drawstrings

Size and packing:

- available sizes: 2.5 cm x 3 cm, 2 cm x 1.5 cm, 3 cm x 2 cm, 3.5 cm x 2.5 cm, 3 cm x 4 cm, 5 cm x 5 cm
- sterile: packed by 1 pc

NOTE: For invasive procedures use the product only in an unwound state.



MD

Devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>