

- made of chlorine-free bleached cotton gauze or medical nonwoven
- all edges folded inside to prevent the release of threads
- small dressing material with high absorption capacity
- different sizes make it possible to adjust the absorption to individual needs
- can be sterilised

Application:

- for draining fluids in surgical and outpatient conditions
- for disinfection of the OP-field

The product is to be used in an unrolled state.

Types:

- sterile and non sterile
- with and without X-ray detectable thread
- 4- and 8-ply

Size and packing:

- sterile: packed by 1 pc, 2 pcs, 5 pcs
- non sterile: packed by 68 pcs, 120 pcs or 200 pcs
- sizes: 1 m x 1 cm, 1 m x 2 cm, 1 m x 5 cm, 2 m x 1 cm, 2 m x 2 cm, 2 m x 5 cm, 2 m x 10 cm, 2 m x 15 cm

NOTE: For invasive procedures use the product only in an unwound state.



MD

Devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>