

- made of chlorine-free bleached cotton gauze or medical nonwoven
- all edges folded inside to prevent the release of threads
- small dressing material with high absorption capacity
- variety of shapes make it possible to adjust the product to individual needs and activities
- different sizes make it possible to adjust the absorption to individual needs
- can be sterilised

#### Application:

- for draining fluids in surgical and outpatient conditions
- for disinfection of the OP-field

#### Types:

- sterile and non sterile
- with and without X-ray detectable thread
- available shapes: balls, beans, whiskers, triangles, thimbles, grains, pepper, for laryngology

#### Size and packing:

- sterile: packed by 3 pcs, 5 pcs, 10 pcs and by 2 x 3 pcs in two chambers of blister packaging (double blister)
- non sterile: packed by 25 pcs, 50 pcs, 100 pcs, 200 pcs, 500 pcs
- wide variety of sizes

NOTE: For invasive procedures use the product only in an unwound state.



**MD**

Devices are certified by the Notified Body TÜV Nord Polska Sp. z o.o. and bear the CE mark with number 2274

>>