

- made of flexible durable fabric that perfectly adjusts to any body shape
- centrally placed absorbent pad is secured with micro-net to protect the dressing from sticking to the wound
- acrylic adhesive does not cause skin irritation
- the length of the plaster enables to adjust the size to individual needs

Application:

- for dressing abrasions, minor cuts and wounds in people with normal skin sensitivity

Types:

- non sterile

Size and packing:

- available sizes: 6 cm x 1 m, 8 cm x 1 m
- non sterile: packed by 1 pc



MD

>>